3X AAV PRODUCTION YIELD USING TFDF® FOR CELL GROWTH INTENSIFICATION AND CLARIFICATION

Rachel Legmann, Earl Pineda, Michael Bransby, Philip Yuen, Rene Gantier; Repligen Corporation, Waltham, MA, USA

Joao Mendes, Ricardo Silva, Cristina Peixoto, Antonio Roldao; iBET, Oeiras, Portugal

Introduction

Gene therapy demand generated by high dosage and large patient populations require significance improvements to viral vector manufacturing. Two integrated perfusion recovery. KrosFlo® TFDF® increased total yield of adeno-associated virus serotype 8 (AAV8) by 3 fold compared to a batch bioreactor process; lysate recovery yield was comparable to depth filtration. This case study illustrates the benefits of intensified and integrated processes for high-productivity and cost-effective viral vector manufacturing to meet the increasing global demand.

AAV Viral Vector Manufacturing Workflow



Viral vector intensification helps meet manufacturing demand

Process step	Batch process
Virus production	 Transfection at ~2E6 cells/mL (~48h after inoculation)
Lysate clarification	2-stage filtration: DF-based processGlass Fiber filter (5 um nominal)
	 0.8/0.2 um (nominal) membrane filter
Inoculation	Seeding of HEK293T cells in BalanCD HEK293 media at ~6E5



Run	Filter	Load (L/m ² @ 0.7bar)	Surface area (m ²) @ 200 L scale	FLUX (LMH)	Clarification step recovery yie	
Batch	Glass fiber (5 um)	23	9	100	70	
	0.8/0.2 um membrane	30	7	100	70	
KrosFlo® TFDF® (Perfusion and clarification)	Integrated TFDF [®] (2 - 5 um)	300	0.6	450	70	
	Multi-layer polypropylene (< 0.5 um)	53	4	100		
	Simplified and smaller for	otprint clarificatio	n process with integrated	KrosFlo® TFD	OF [®] filter	

Conclusion

• Krosflo[®] TFDF[®] provides one solution for upstream viral vector intensification and clarification • High-productivity and cost-effective manufacturing processes help meet the global demand for viral vectors

© 2022 Repligen Corporation. All rights reserved. The trademarks mentioned herein are the property of Repligen Corporation and/or its affiliate(s) or their respective owners.



Intensified process (XCell ATF[®] or KrosFlo[®] TFDF[®])

- Transfection at ~5E6 cells/mL (perfusion started 48 hours after inoculation at 1 vvd)
- Integrated first stage clarification with TFDF[®] (2 5 um nominal)

cells/mL in glass vessel bioreactor (2 L working volume)



- > 3-fold increased AAV8 production yield
- Improved process economics

• Simplified process enhances cell specific productivity and growth leading to higher vector yield

25 batches, the cost savings per year can

be almost \$24M



Patients /year	vg*/year	Cell culture mode	Cell cult. vol. /year (KL)	2 KL batches/ year	# of 2KL bioreactors @ 25 batches/year	Consumable cost*/year (M\$)	Cost* savings/ year (M\$)
5000	6 x 10 ¹⁸	Batch	148	74	3	72	-
5000	8 x 10 ¹⁹	Perfusion	50	25	1	48	24
		Batch	980	490	20	482	-
		Perfusion	330	166	9	320	162